

STATIONARY OXYGEN CONCENTRATOR



* Humidifier bottle and canula sold separately

ITEM #	DESCRIPTION
LM5A	5 LPM Stationary Oxygen Concentrator

FEATURES

- Compact footprint with top handle for easy transportation and storage.
- Reliable and safe. High efficiency.
- Easy-to-read digital hour meter display.
- Safety Auto shutdown for power surges, high and low pressure, high temperature, and low oxygen concentration.
- Heat protection to ensure safety of the compressor and concentrator.
- Unit comes standard with oxygen purity and low flow rate indicators.

SPECIFICATIONS

Size	13" W x 10.3" D x 21.3" H
Weight	35.5 lbs
Power Consumption	360W (average)
Alarms	Power failure, sensor failure, <73% ox lasting 10 min.
O ₂ Purity	93%± 3%*
Outlet Pressure	5 PSI ± .73 PSI
Pressure Relief Mechanism	36 PSI ± 7.3 PSI
Maximum Recommended Flow	5L
Liter flow range	0.5~5L/min ,93%±3% (after running for 30 minutes)
Power Supply:	AC120V, 60Hz; Current:3.5A; Power: 450VA
Use	Designed for 24/7 use

V3: 12/16/20