

# POC Comparison Chart



MANUFACTURER	Rhythm Healthcare	Inogen One <sup>®</sup>	Phillips
MODEL	P2	G3	SimplyGo Mini
Height	6.43"	8.75"	9.4"
Width	3.4"	3"	3.6"
Depth	8.7"	7.25"	8.3"
Weight (inc. battery)	4.37 lbs	4.7 lbs	5 lbs
Hours of Use (inc. 1 battery)	Up to 5 hours on setting 1	Up to 4.5 hours on setting of 2	Up to 4.5 hours on setting of 2
Patient Changeable Sieve	No	Yes	No
Pulse Settings	1-5	1-5	1-5
Warranty	3 Years	3 Years	3 Years