


PERFECT FOR PATIENTS
PERFECT FOR CAREGIVERS



CREATING INNOVATION AND DIFFERENCE

Rhythm·Turn

Rhythm·Turn offers continual lateral rotation therapy to help mobilize secretions and other bodily fluids during pulmonary treatment related to immobility.



Various top covers are designed and available to satisfy your individual requirement in terms of therapy, prevention and patient comfort.

Master Control Unit

- Dimension: (L)17.7"x (W)6.8"x (H)10.8"
- Weight: 13.8 lbs
- Rated Voltage: 110-120 volts
- Max Current: 5A
- Rated Frequency: 60Hz
- Cycle Time: 3-95 mins
- 3 in 1 alternation
- Min/Max Pressure: 10-61 mmHg +/- 3-4 mmHg
- Protection Type: Class I, Type BF
- Power outage visual and audio alarm
- Static Function
- Auto Firm for fast mattress inflation
- Keypad Lock-out Function

Mattress Replacement

- Model No: Rhythm·Turn
- Dimension: (L)80"x (W)36"x (H)10"
- Weight: 37.4 lbs
- Weight Capacity: 600 lbs
- Cell Material: Blue Nylon with PU Coating
- Cell Number: 18 cells
- Cover Material: Blue Nylon with PU coating finish
- Base Material: Black Polyester with PVC Coating
- CPR Quick Release System

Toll Free:
877-367-5507

LA, CA OFFICE
TEL: 626-275-4063
TAMPA, FL OFFICE
TEL: 813-854-5464



Rhythm·Turn



Rhythm·Turn

provides you, the carer, with peace of mind
and your patient with world class care & comfort

Perfect for patients and caregivers

- High vapor permeable and quilted nylon cover

- Low Air Loss Therapy provides extra comfort for moisture and temperature control

- The Rhythm-Turn replacement mattress can accommodate patient weights up to 600 lbs

- 20°C and 40°C lateral rotation with side perimeter and static therapy



Low slip base cover is a plus to keep the mattress in place on the bed



Effective micro low air loss for moisture management



Nylon woven fabric with PU coating finish

Rhythm-Turn is your partner in protection and treatment of pressure ulcer



CPR quick release system deflates the system within 20 seconds.



Cable management keeps electrical cables off the floor



Hanging bracket fits on all type of beds

